



whealthTM Spiritus-5

Robust * Responsive * Reliable

Spiritus-5

1-5lit Oxygen Concentrator



* whealth Spiritus-5 is an extremely reliable device with benchmark quality & purity of continuous oxygen supply

* Meticulously designed and made by understanding needs of the customer; Spiritus-5 is easy to use. One button start makes it user friendly for elder patients

* Ideal for patient with Lung Disorders like ILD, COPD, PF and having low oxygen saturation



whealthTM

Specifications:



Model	Spiritus-5
Oxygen Flow	0.5 - 5 Lit/ Min
Oxygen Purity	93% (+/- 3%)
Outlet Pressure	50 - 70 kPa
Flow meter	Single
Therapy Hour Indicator	Yes
Alarms	Power Failure High / Low Pressure Low Oxygen Purity
Decibel Level	40 dB
Net Weight	13.5 Kg
Portable	vide castor wheels
Power Supply	220V / 50 Hz 110V / 60Hz
Power Consumption	300 W
Power Cable	India
Electric Safety	Protection Class II, Type B
Warranty	3 years - Serviceable

Attractive 3 Years
Warranty

Low Maintenance

Spectacular Oxygen
Purity



whealth[™]

Marketed By:

**Evron Biocraft
Pvt Ltd**

B-11, Chand Soc,
Ground Floor,
Alfred Creado Road, Juhu,
Mumbai-49
+91-22-26255502/ 08



www.evron.in

care@evron.in

whealth[™]

Service Centre:

B-11, Chand Society, Ground Floor, Alfred Creado Road, Near ISKCON temple, Next To Vaishali Hotel, Juhu, Mumbai-400049.

Cust Care No: 022-26255508 . Email : service@evron.in

Evron Biocraft Pvt Ltd reserves the right to make changes in specification and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: Device sale is restricted to sale by or on the order of a physician