




## MULTI-FUNCTIONAL OXYGEN CONCENTRATOR

Focus On Medicalcare O2 Solutions!

All For Stability



### Recommendation

- The new EVOX-5 Version is a SmartTouch multifunctional oxygen concentrator
- It has been specially developed to provide a dependable, convenient source of supplemental oxygen for therapeutic needs.
- Includes Nebulizer function, Individual nebulizer outlet
- Large LED Display
- Oxygen Purity Display

Suitable for children and adults, EVOX-5 offers concentration up to 96% and flow rates reaching 5LPM.



Technical specifications

System Map	Function	EVOX-5S	EVOX-5
Oxygen Generating System	Oxygen Flow	0.5~5L/min	0.5~5L/min
	Oxygen Concentration	93%+3%	93%+3%
	Oxygen Output Pressure	0.04-0.07Mpa	0.04-0.07Mpa
	Oxygen Number Display	YES	NO

Nebulize System	Nebulize Mode	YES	-
	Individual Nebulizer Outlet	YES	-
	Display Method	On LED screen	-
	Operating Way	By touch panel button	-
	Accessories	Nebulize nozzle and mask	-
	Nebulizing particle Size	3.6um+25%	-

Screen Display System	Display Material	All by LED lights
	What Can Be Displayed ?	“Oxygen” light on
		“Nebulize” light on
		“Oxygen Percentage” light on
		“Each Runing Time” light on
		“Accumulated Time” light on
		“Timer Setting Number”

Operating System	Operating Panel	Smart Touch
	OXYGEN	Start to produce Oxygen
	NEBULIZE	Start to run Nebulizer
	TIMER“+”and “-”	Seting from 0 to 3 hours (15mins increase or decrease)

Safety System	Low Oxygen Alarm	Yes
	Oxygen Pressure Protection Alarm	Yes
	Pressure Faulty Alarm	Yes
	Power off Alarm	Yes
	Compressor Faulty Alarm	Yes
	Overheated Protection	Yes

Electrical System	Power Consumption	350 W
	Auto Power off by timer	YES
	Screen Sleep Mode	YES
	Noise Level	<43dB(A)

Packing	Machine Body Size	380x270x580mm
	Package Size	460x350x680mm
	Net Weight	15.6kg
	Gross Weight	19kg
	Body Material	ABS, Glaze

Optional Functions	Remote Control	Optional
	Voice Guide System	Optional