

INVACARE®
PERFECTO₂™
OXYGEN CONCENTRATOR



The **Invacare® Perfecto₂™ Oxygen Concentrator** is the smallest, lightest, quietest and most efficient 5-liter concentrator ever produced by Invacare. The high performing Perfecto₂ unit incorporates all the best features of the market-leading Invacare® Platinum™ XL concentrator, while significantly improving upon the key specifications of the Platinum XL Concentrator and offering a contemporary, patient-preferred design and color. Available with or without the **Invacare® SensO₂® Oxygen Sensor**, the Perfecto₂ Concentrator represents a quantum leap forward in the evolution of oxygen concentrators.

Features

- 75% quieter than Invacare® Platinum™ XL Concentrator
- 25% more energy efficient than Invacare Platinum XL Concentrator
- 17% lighter than Invacare Platinum XL Concentrator
- 33% smaller than Invacare Platinum XL Concentrator
- Contemporary design and color
- Top handle for ease of mobility
- Pressure-based system is superior to timing-based systems

- Self-diagnostic electronics minimize troubleshooting time
- Servicing parts compatibility: Compatible with many Platinum concentrator parts
- Alarms and safety systems: Battery-free power-loss alarm; Invacare® Sieve-GARD™ Monitor protects sieve from degradation; Invacare® SensO₂® Oxygen Monitor; Compressor 35 PSI pressure-relief valve, Low flow alarm
- Simplified preventative maintenance schedule: Replace preventative maintenance components as needed

HCPCS Code: E1390



Invacare® HomeFill® Oxygen System connection
 HomeFill compressor and Ready-Rack have been updated to coordinate with the Perfecto₂ Concentrator. Use with Invacare cylinders
 Easy to maintain filter



Invacare Corporation
 www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2012
 Invacare Corporation.
 All rights reserved.
 Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
 Form No. 07-068
 Rev. 01/12
 111128

Model #	Dimensions	Weight	Shipping Weight	Quantity	Oxygen Flow	Flow Range	Maximum Outlet Pressure	Pressure Relief Mechanism
IRC5P/ IRC5PO2	23" H 13" W 11.5" D	43 lb.	47 lb.	1 Per Case 27 Per Pallet	5 LPM	0.5-5 LPM	5+/- 0.5 PSI	35 PSI +/- 3.5 PSI
	O₂	Filter	Electrical Requirements	Rated Current Input	Average Power Consumption	Sound Level	Altitude	Limited Warranty
	95.6%-87% At All Flow Rates	Cabinet Compressor Inlet Outlet HEPA	120 VAC +/- 10%, 60 Hz	4.0 A	280 W @ 3 LPM	43 db.	8,000 ft. 6,000 ft. w/ HomeFill System	3 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

