



# whealth<sup>TM</sup> Spiritus-10

Trust in the Best

Spiritus-10

10 lit Oxygen Concentrator



\* whealth Spiritus-10 is an extremely reliable device with benchmark quality & purity of continuous oxygen supply

\* Meticulously designed and made by understanding needs of the customer; Spiritus-10 is easy to use. One button start makes it user friendly for elder patients

\* Ideal for patient with Lung Disorders like ILD, COPD, PF and having low oxygen saturation due to Covid-19 like symptoms



whealth<sup>TM</sup>


## Specifications:



Large LCD



Timer



Oxygen



Pressure



Atomizer



Flow

Model	Spiritus - 10
Oxygen Flow	2 - 10 Lit/ Min
Oxygen Purity	93% (+/- 3%)
Outlet Pressure	50 - 70 kPa
Outlet for Oxygen	Single (Double outlet is optional)
Atomizer	Optional; on select models only
Flow meter	Single
Alarms	Power Failure
	High / Low Pressure
	Low Oxygen Purity
Decibel Level	50 dB
Net Weight	28 Kg
Portable	vide castor wheels
Power Supply	220V / 50 Hz
Power Consumption	720 W
Power Cable	India
Electric Safety	Protection Class II, Type B
Warranty	1 years - Serviceable

Attractive 1 Year  
Warranty

Low Maintenance

Spectacular Oxygen  
Purity



whealth™

Marketed By:

**Evron Biocraft  
Pvt Ltd**

B-11, Chand Soc,  
Ground Floor,  
Alfred Creado Road, Juhu,  
Mumbai-49  
+91-22-26255502/ 08



[www.evron.in](http://www.evron.in)

[care@evron.in](mailto:care@evron.in)

whealth™

**Service Centre:**

B-11, Chand Society, Ground Floor, Alfred Creado Road, Near ISKCON temple, Next To Vaishali Hotel, Juhu, Mumbai-400049.

**Cust Care No: 022-26255508 . Email : [service@evron.in](mailto:service@evron.in)**

Evron Biocraft Pvt Ltd reserves the right to make changes in specification and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

**CAUTION:** Device sale is restricted to sale by or on the order of a physician