

Esleh Grand



Features:

- Foldable Electric
- Flip up Armrest
- Adjustable Footrest
- Removable battery For easy Charging
- Wireless remote access (Optional)
- Self propelling options

• Technical specifications:

- Battery 24v V 24 AH (lead acid double Battery) (lithium ion battery optional)
- Electric brake
- Motor : DC 250 W 2 Pcs
- Frame : Steel (3 layer Powder coating)
- Speed : 6 km/h
- Driving Range : 15-20 Km (in radius)
- Rear wheel : 22 inch
- Front wheel : 10 inch
- Seat Height : 18 inches (floor)
- Seat width : 18 inches
- Seat depth : 18 inches
- Overall width : 24 inches
- Backrest height : 18 inches
- Overall height : 48 inch
- Total wheelchair weight : 39 kgs
- Battery weight : 7 kgs
- Load bearing : standard 120 kgs, max 150 kgs