

Esleh Master



Features:

Reclining backrest

- Easy Folding
- Forest position can be Adjusted by remote
- Wireless remote access
- Control by attendant
- Flip up armrest
- Controller can be operate from both side
- Head Support
-

Technical specifications:

- Battery 24V 24 AH (lead acid double Battery)
- Electromagnetic brake
- Motor : DC 250 W 2 Pcs
- Frame : Steel (3 layer Powder coating)
- Speed : Speed: 0-6km/h, five gear can be adjusted
- Driving Range : 15-20 Km (in radius)
- Controller : 360 degree rocker controller
- Front wheel : 10 inch
- Rear wheel : 16 inch
- Seat Height : 19 inch
- Seat width : 18 inch
- Seat depth : 18 inch
- Overall width : 24 inch
- Backrest height : 18 inch
- Overall height : 48 inch
- Total wheelchair weight : 49 kgs
- Battery weight : 7 kgs
- Load bearing : standard 120 kgs, max 150 kgs