

Esleh Super Plus



Features:

- Light Weight
- Light Indicators
- Easy in Air Traveling
- Quick Folding
- Adjustable hand-rest
- Wireless Remote Access
- Control by Attendant
- Reclining Backrest
- Controller can be operate from both side
- Comfortable Head support
- Anti tilting wheels, safe and reliable
- Safety belts
- Storage bag
- Frame : Aluminium
- Alloy : Aluminium
- Breathable cushion
- Large area foot pedal to relieve muscle pressure
- Off road wheel
- Removable battery
- Smart controller for simple operation
- Load bearing : standard 120 kgs, max 150 kgs

Specifications :

- Battery 24 V 20 AH (Lithium ion Battery)
- Electromagnetic Brake
- Motor: DC 250 W 2 Pcs
- Powerful Brushless motor
- Automatic parking
- Speed: 0-6km/h, five gear can be adjusted
- Driving Range : 20-25 kms (In radius)
- 360 degree rocker controller
- Shock absorption to reduce vibration over harsh terrain
- Endurance +,- 25km
- Rear wheel : 12 inch
- Front wheel : 8 inch
- Seat Height : 18 inch (floor)
- Seat width : 18 inch
- Seat depth : 18 inch
- Overall width : 24 inch
- Backrest height : 18 inch
- Overall height : 37 inch
- Battery weight : 2.2 kgs
- Total wheelchair weight : 26 kgs

Wireless Remote (Optional with Each Model except Stair chair)