

Esleh Travel



Features:

- Adjustable hand-rest
- Controller can be operate from both side
- Extra control pannel
- Joystick with 2 inch large display
- Anti tilting wheels, safe and reliable
- Safety belts
- Storage bag
- Frame : Aluminium
- Alloy : Aluminium
- Battery weight : 2.2 kgs
- Breathable cushion
- Large area foot pedal to relieve muscle pressure
- Off road wheel
- Removable battery
- Smart controller for simple operation
- Load bearing : standard 120 kgs, max 150 kgs

Specifications :

- Battery 24 V 20 AH (Lithium ion Battery)
- Electromagnetic Brake
- Motor: DC 250 W 2 Pcs
- Frame : Aluminium
- Speed 6 km/h
- Driving Range : 20-25 kms (In radius)
- Rear wheel : 9.5 inch
- Front wheel : 7 inch
- Seat Height : 18 inch (floor)
- Seat width : 16 inch
- Seat depth : 16 inch
- Overall width : 21.5 inch
- Backrest height : 18 inch
- Total wheelchair weight : 19 kgs
- Battery weight : 2 kgs
- Load bearing : standard 120 kgs, max 150 kgs