


# INOGEN ONE<sup>®</sup> G5<sup>™</sup> System



Inogen One  
AC Power Supply\*  
Item #BA-501



Inogen One G5  
DC Power Cable\*  
Item #BA-306



Inogen One G5  
Lithium Ion Battery\*  
Up to 6.5 hours run time  
Item #BA-500



Inogen One G5  
Carry Bag\*  
Item #CA-500



Inogen One G5  
Backpack\*\*  
Item #CA-550



Inogen One<sup>®</sup> G5  
Extended Life Lithium Ion Battery\*\*  
Up to 13 hours run time  
Item #BA-516

## The Inogen One G5 System (IS-500):

- Compact and Lightweight
- Flow Settings 1-6
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use
- Easy-to-read LCD display
- Inogen Connect App

\*Included with system

\*\*Accessory sold separately

## Oxygen. Anytime. Anywhere<sup>®</sup>

Once again, Inogen innovation sets a new standard in oxygen delivery.

The Inogen One G5 offers the most oxygen per pound for a portable oxygen concentrator on the market today. With flow settings from 1-6, the Inogen One G5 is designed to dramatically increase independence for most users 24/7. The Inogen One G5's user interface has been designed with an easy-to-read LCD display and simple control functions. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere.

The Inogen One G5 is your one solution, your single solution, for oxygen at home or away.

|             |                                                                                                                                                                                     |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Weight      | 4.7 lbs. (includes single battery)                                                                                                                                                  |
| Size        | Length: 7.19"<br>Width: 3.26"<br>Height: 8.15" (includes single battery)                                                                                                            |
| Oxygen Flow | Pulse dose delivery system<br>Six flow settings 1, 2, 3, 4, 5, 6                                                                                                                    |
| Power       | AC Power Supply 100-240V, 50-60Hz<br>(auto sensing to allow worldwide use)<br>DC Power Cable: for mobile use in car                                                                 |
| Battery     | Duration (single battery): Up to 6.5 hours<br>Recharge up to 3 hours with AC or DC power<br>Duration (double battery): Up to 13 hours<br>Recharge up to 6 hours with AC or DC power |
| Noise Level | 38 dBA*                                                                                                                                                                             |
| Warranty    | 3 year or lifetime limited warranty**                                                                                                                                               |
| Operation   | Simple control functions and easy-to-read LCD display<br>Inogen Connect App                                                                                                         |
| Use         | Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft                                                                              |

\*At flow setting 2

\*\*See Inogen One<sup>®</sup> Limited Warranty Statement for complete warranty information

Manufactured by:

**inogen**  
326 Bollay Drive  
Goleta, California 93117  
[www.inogen.com](http://www.inogen.com)



## Marketed by

Office No. 506, 5th Floor, Office Tower – B, Elante,  
Plot no. 178-178A, Ind Area, Phase I, Chandigarh

Call : +91 7341103371  
[www.sanraishop.com](http://www.sanraishop.com)