




BPAP PLUS SERIES



Shenyang RMS Medical Tech Co. Ltd

Picture for reference only is subject to material object

Product Specification

Model	BPAP20 PLUS	BPAP25 PLUS	BPAP30 PLUS
Modes	CPAP,S,T,S/T	CPAP,S,T,S/T	CPAP,S,T,S/T
IPAP range	4-20 CmH2O	4-25 CmH2O	4-30 CmH2O
EPAP range	4-20 CmH2O	4-20 CmH2O	4-20 CmH2O
CPAP range	4-20 CmH2O	4-20 CmH2O	4-20 CmH2O
I-trigger	Auto,1-6	Auto,1-6	Auto,1-6
E-trigger	Auto,1-6	Auto,1-6	Auto,1-6
IPAP rise time	Auto,0.1-0.6s	Auto,0.1-0.6s	Auto,0.1-0.6s
Breath rate	2-40 bpm	2-40 bpm	2-40 bpm
Noise Level (dB)	<28dB	<28dB	<28dB
Ramp time(minute)	0-60	0-60	0-60
Display data	Vt,Vleak,Rate,I:E	IPAP, EPAP, Rate(bpm), I:E, Vt, Vleak, Mvent, Trigger I/E	IPAP, EPAP, Rate(bpm), I:E, Vt, Vleak, Mvent, Trigger I/E
Alarms	Tube alarm, Apnea alarm & power failure	Tube Alarm, Apnea Alarm, High/ low rate Alarms & power failure	Tube Alarm, Apnea Alarm, High/ low rate Alarms Mvent alarm & power failure
Dimensions	18*27*11cm	18*27*11cm	18*27*11cm
Weight	1.7 Kg	1.7 Kg	1.7 Kg
Screen size	2.4 inch	4.3 inch	4.3 inch
Heated humidifier	Integrated	Integrated	Integrated
Memory	Micro SD Card - 8 GB	Micro SD Card - 8 GB	Micro SD Card - 8 GB
Warranty	2 years	2 years	2 years

Electrical Characteristics

Electrical supply	100-240V, 50-60 Hz
Power Consumption	≤ 80W with Humidifier
93/42/EC Classification	Class II A

Environmental requirements

Temperature	5°C-35°C
Humidity	30-95%, no condensation
Storage	-20-60-°C
Altitude	From Sea level to 2700 mtrs



Shenyang RMS Medical Tech Co. Ltd.

India Office :


Mobile: +91-9717986888, 9717228868

E-mail : pankaj@rms-care.com

Head Office: No. 21, Gaoke Road,

Hunnan District, Shenyang. Zip Code - 110169

www.rms-care.com/en



Mode->S/T		
Mode	S/T	
IPAP	4.0	CmH2O
EPAP	0.40	CmH2O
Rise Time	0.1	Sec
Rate	10.0	bpm
Insp Time	1.5	Sec
I-trigger	6	